

Note that this data sheet is not lot-specific. Please consult the vial label and the certificate of analysis for information on specific lots.

B-TeZ IAC T2/HT2 3 ml

Catalogue Number: BTTE316005

Package Size: 10 Columns

1. Description

Immuno affinity column based on monoclonal antibodies bonded to macroporous pearl cellulose beads.

2. Specificity

T2-Toxin and HT2

3. Form

Columns filled with gel in phosphate buffered saline (50 mM K-phosphate, pH 7.5, 150 mM NaCl) containing 0.05% sodium azide and 0.01% 2-methyl-4-isothiazolin-3-one embedded within two porous frits.

4. Gel Bed Volume:

0.3 - 0.5 ml

5. Column Working Range:

1⁽¹⁾ – 500ng T2/HT2 ⁽¹⁾Lower limit only depends on the sensitivity of the detection system (LC-FLD, LC-MS); the given value represents sensitivity of a common LC-FLD system after pre-column derivatization using anthrolyl nitrile/4-DMAP.

6. Recoveries:

85 – 110% Of each T2-species, given under 2. within the column working range

7. Organic Solvent Tolerance:

Up to 15% methanol or 10% acetonitrile in loading solution

8. pH Tolerance:

6.5 – 8.5

9. Flow Rate:

1 – 3ml/min

Note that this data sheet is not lot-specific. Please consult the vial label and the certificate of analysis for information on specific lots.

10. Elution:

1.5- 3 ml Methanol

11. Column Capacity:

3000ng T2-Toxin

12. Cross Reactivity⁽²⁾:

T2-Toxin (100%)

HT2 (118%)

Deoxynivalenol (<1%)

⁽²⁾Relative recovery rates when a semi-saturative quantity of T2-Toxin/HT2 and DON (1500ng) in mixture is applied on the column.

13. Application:

Clean up of ailment (and feeding stuffs) extracts for reliable identification and accurate quantification of T2-toxin.

14. Storage

The performance of the immuno affinity column is guaranteed until the expiry date given on the label or certificate if stored at 4°C. Avoid freezing of the product.

15. Additional Information

Columns should be free of air bubbles inside the gel to achieve optimal results regarding flow and recovery rates.

16. Product Instruction

This product is shipped along with an instruction where detailed information is given how to use B-TeZ immuno affinity columns at its best.