

Note that this data sheet is not lot-specific. Please consult the vial label and the certificate of analysis for information on specific lots.

Anti-TA Antibody (PAS Antibody)

Catalogue number: 12 700 010 Package size: 1.0 mg

Overview

Antigen	TA
Type	Primary
Clonality	Monoclonal
Clone	#BCAB1400
Host	Mouse
Application	ELISA, Lateral Flow Assay
Quantity	1.0 mg
Order number	12 700 010

Additional Information

Format	Lyophilized
Handling Advice	Reconstitute antibody by adding 1 ml deionized water.
Isotype	IgG1
Specificity	100% TA
Molecular Weight	150000 Da

Application Details

Application Notes	For capturing of TA labelled antibodies or other proteins in ELISA and Lateral Flow Assay. Can be used for coating. Each laboratory should determine an optimum working concentration for use in its particular application. Other applications have not been tested but they should not necessarily be excluded..
Purity	>95% pure. Purity is tested by electrophoresis.
Purification	Protein A purified
Buffer	Neutralized citrate buffer, pH=7.4
Preservative	0.05 % sodium azide. Warning: This product contains sodium azide, which has been classified as harmful in the concentration range of 0.1–1.0 %.
Storage	Long-term storage: -20° C, Short-term storage: 2-8° C.
Restrictions	For research use only.
Related products:	TAPAS Labelling Kit BTTAP-001 or -003 for 1 / 3 Labels