

D A T A S H E E T

P A S A n t i b o d y P o l y - H R P

PAS Antibody, Poly-HRP labeled

Order number:	20 700 001
Target:	PAS antibody (monoclonal anti-TA antibody, PAS Antibody binds highly efficient the TA-Tag)
Label:	Poly-HRP Poly-Peroxidase labels using highly activated Poly-HRP (Horseradish Peroxidase). Poly-Peroxidase labels are characterized by high measurement signal amplification.
Purification	Size-exclusion chromatography
Conjugate Quantity:	200 µg
Format:	HRP stabilizer solution (contains preservative: 0.02% methylisothiazolone, 0.02% bromonitrodioxane and 0.02% other active isothiazolones
Application:	ELISA / IHC PAS Antibody Poly-HRP labelled binds highly efficient the TA-Tag and provides a very strong signal. The conjugate supports the TA-PAS Labelling System Recommended ELISA concentration: approx. 25 ng/ml Each laboratory should determine the optimum working dilution for use in its particular application.
Storage:	2-8°C, Keep dark.
Shelf Life:	>6 months
Restrictions:	For research use only.